



Nuclear Medicine


Sahar Rezaei, PhD

Assistant Professor
Tabriz University of Medical Sciences

PhD in Medical physics, Department of Medical physics and Biomedical Engineering, Tehran University Of Medical Sciences


Medical Imaging

- Show and record the structural or functional state of the body.
- that cannot be observed by visual inspection or physical examination




Medical Imaging -Image Matrix

- Image generation begins with an empty matrix.
- Image matrix is filled by processed acquired data.



Medical Imaging -Image Matrix

We need to obtain each pixel value in imaging process




Medical Imaging

Imaging is the creation of a 2D, 3D or 4D representation of a physical phenomena by measurement of energy

- Emitted
- Reflected
- Transmitted

from the object.

- Structural/Anatomical Imaging
- Functional Imaging
- Hybrid Imaging(mostly Anatomical+Functional)



Structural Imaging

- Radiology
- Fluoroscopy / Angiography
- Computed Tomography (CT)
- Mammography
- MRI
- Ultrasound



Functional Imaging

- Gamma Camera
- SPECT
- o PET
- o fMRI
- Molecular Imaging



Hybrid Imaging


- SPECT/CT
- PET/CT
- PET/MRI
- MRI/Mammo

. . .




What is Nuclear Medicine?

Nuclear medicine is a branch of medical imaging



How are images made


- Patient receives radioactive tracer tagged to a chemical
 - via injection, inhalation, or ingestion
 - painless process, cannot feel the painless process,
 - Small radiation dose equivalent to x-ray





How are images made

- Radioactive chemical travels to organ of interest
 - type of chemical determines into which organ(s) the radioactivity distributes




Whole body bone scan showing increased uptake in the lumbar Spine, maxillary sinuses, hemipelvis, right femoral neck, ribs

How are images made

Patient imaged with camera


- as radiation decays, it emits energy that exits the body
- o camera absorbs this energy and creates a picture




Nuclear Medicine Instrumentation

- Gamma Scintillation Camera
- Single-Photon Emission Computed Tomography (SPECT)
 - SPECT/CT
- Positron Emission Tomography (PET)
 - PET/CT

Gamma Scintillation Camera




Gamma Scintillation Camera



Nuclear Medicine Instrumentation

- Gamma Scintillation Camera
- Single-Photon Emission Computed Tomography (SPECT)
 - SPECT/CT
- Positron Emission Tomography (PET)
 - PET/CT


SPECT












SPECT




Dual Head Gamma Camera in Cardiac Configuration

Characteristics of Commonly Used Radionuclides

	Symbol	Physical Half-Life	Approximate Energy
Photon-Emitting Radionuclides for Imaging	-		Gamma (keV)
Technetium-99m	^{99m} Tc	6 h	140
Molybdenum-99	⁹⁹ Mo	67 h	181, 740, 780
lodine-123	123	13.2 h	159
lodine-131	¹³¹	8.0 days	364
Xenon-133	¹³³ Xe	5.3 days	81
Gallium-67	⁶⁷ Ga	78.3 h	93, 184, 296, 388
Indium-111	¹¹¹ In	67 h	173, 247
Indium-113m	^{113m} In	1.7 h	392
Thallium-201	²⁰¹ TI	73.1 h	69, 81 (x-rays from mercury daughter)
Krypton-81m	^{81m} Kr	13 s	191


SPECT/CT




SPECT/CT

Benefits of SPECT/CT

Functional mapping






SPECT/CT

Benefits of SPECT/CT

- Functional mapping
- Quantitative SPECT

SPECT & SPECT/CT

- Cardiac Imaging
- Inflammation Imaging
- Brain Imaging
- Radiation Planning



A 60-year-old female patient with cardiac sarcoidosis.

Nuclear Medicine Instrumentation

- Gamma Scintillation Camera
- Single-Photon Emission Computed Tomography (SPECT)
 - SPECT/CT
- Positron Emission Tomography (PET)
 - PET/CT



Half life

20 minutes

10 minutes

2 minutes

8 minutes

10 minutes

12.7 hours

68.1 hours

1.3 minutes

4.2 days

110 minutes

Positron energy (MeV)

0.385

0.492

0.735

0.250

1.216

1.315

0.278

0.836, 0.352

1.523, 1.157

1.691, 7.228, 1.509, 1.376

P	E.	I	Ra	di	lon	ıu	cli	de	3. (Ç

Isotope

C - 11

N-13

O-15

F-18

K-38


Cu-62

Cu-64

Ga-68

Rb-82


I-124






PET radiopharmaceuticals for brain imaging


Compound	Application
O-15 H ₂ O	Blood flow
O-15 O ₂	Oxygen metabolism / flow
O-15 or C-11 Carboxy haemoglobin	Blood volume
C-11 Methionine	Amino acid metabolism
C-11 Ephedrine	Adrenergic terminals
C-11 Carfentanil	Opiate receptor activity
C-11 Flunitrazepam	Benzodiazepine receptor activity
C-11 Methylspiperone	Dopamine receptor activity
C-11 Scopolamine	Muscarinic cholinergic receptor activity
F-18 Fluorodeoxyglucose (FDG)	Glucose metabolism
F-18 Fluoro-dopa	Presynaptic dopamine system
F-18 Fluorothymidine (FLT)	DNA synthesis





PET/CT





Hybrid PET/CT scanner.

















Clinical Applications of PET and PET-CT




FDG PET/CT for staging a case of NHL




FDG uptake at the periphery of an abdominal mass lesion with central cold area of necrosis




Brain FDG PET/CT scan



FDG PET/CT Scan for solitary pulmonary nodule



FDG PET/CT for therapeutic assessment in a case of metastatic breast carcinoma



Cardiac PET/CT Scan

End of this Part